

CoviCare EXECUTION GUIDE

CONTENTS

- **Introduction**
- **Prerequisites**
- **Hardware**
- **Hardware setup**
- **IBM account**
- **How to use CoviCare ?**
- **Conclusion**

Introduction

The covid pandemic has completely changed our outlook towards everything. Post covid, there is going to be massive changes in the work culture of people. Bringing a change in the way people work is the absolute need of the hour. Companies cannot afford to overlook the health conditions of their employees. They also have to ensure that their employees practice social distancing. Even if one employee in the company is infected, it can easily spread to others. This implies that workplaces can turn out to be potential hotspots if things are left unchecked. Hence, **Social distancing and constant health monitoring is required in workplaces where thousands work.**

Our idea is to implement a lost cost solution that can bring about social distancing and can track the health of the employees at the same time but not at the cost of their valuable time. It ensures a normal working environment with a minimal change.

- Node MCU
- IBM Cloud
- Ultrasonic Sensor HC-SR04
- Pulse Sensor
- MLX90614
- Nodered
- GPS Module NEO-6MV2
- RFID Tag
- Buzzer

Prerequisites


Software

- Arduino IDE




How to install Arduino IDE?

Step 1-Download the most suitable version from the Arduino official website. Unzip the file.




Step 2: Power up your Nodemcu. Inside the unzipped folder you can find the application icon with infinity label. Double-click to start the IDE.



ARDUINO 1.8.12 is installed.


- Library for NodeMCU

Link: "<https://github.com/esp8266/Arduino>"




Library for NodeMCU is installed from the given link as shown above. This Library has to be used for implementing the features of the project.

To include the library:- Sketch->include library->Click on the required one or Add .ZIP library.



- Tinygps+ (for NEO 6M module) Library



Install the Tinygps+ library from the link

https://github.com/codegardenllc/tiny_gps_plus to include the library used to implement the gps tracking feature which would help determine the location of the worker.

Arduino Environment

In tools menu set the following:

- Board : NodeMCU 1.0 (ESP-12E Module)
- Upload Speed : 115200
- CPU Frequency : 160 MHz
- Flash Size : 4M (1M SPIFFS)
- Debug Port : Disabled
- Debug Level : None
- lwIP Variant : v2 Lower Memory
- vTables: Flash
- Exceptions: Disabled
- Erase Flash : Only Sketch
- SSL Support : Basic SSL ciphers (lower ROM use)
- Port : *Connect the ESP8266 to your laptop using a MicroUSB cable and then select your port, depending on OS*

Once the environment is set up, choose the code file, compile it and then flash it to the board. We will also need a local WiFi network to connect nodemcu.

Also add the following lines to the arduino code.

```
// Add WiFi connection information
```

```
char ssid[] = " "; // your network SSID (name)
```

```
char pass[] = " "; // your network password
```

Hardware


The hardware required is as follows:-

1. Nodemcu
2. RFID Tag
3. Ultrasonic sensor HC-SR04
4. Pulse Sensor
5. MLX90614
6. GPS


7. Push button
8. Jumper wires

HARDWARE SETUP

Circuit of module with user:




Scanner:



IBM ACCOUNT


Create IBM Account:

Go to [IBM Platform](#) to create an IBM account.



To create Node-Red Application

- 1.From Dashboard select catalog.
- 2.Select Software
- 3.Select Web and Application
- 4.Select Node-Red App




5. Keep all the values to default and continue with the App development.

- Wait until the Cloudant service has been deployed
- Press the Create Service + button, then select the Internet of Things section and press the Next button.
- Select the Internet of Things Platform then press the Next button
- Press Create to add the Internet of Things Platform to your application.
- Also generate the IBM cloud API key.

6. Once you select on click **Create app**, then proceed to click on **Deploy App**.

7. After deployment wait for the status to turn into **Success**.

8. App URL will appear after successful deployment. Click on the App URL



9. Add connections to the associated services clicking Connections --> Add Connections.

10. Click on Visit App URL in the Overview page.

The screenshot shows the IBM Cloud console for a Node RED KKHIS instance. The top navigation bar includes the IBM Cloud logo, a search bar, and links to Catalog, Docs, Support, and Manage. The main content area is divided into several sections:

- Overview:** Shows the instance is **Running** with a **Visit App URL** link. A red box highlights the **Running** status and the **Visit App URL** link.
- Instances:** Displays the health status as **100%** (1/1 instance(s) are running) and a slider for MB memory per instance (0 to 256).
- Runtime:** Shows the SDK for Node.js™ and a donut chart indicating **256** Total MB allocation, with 0 MB still available.
- Runtime cost:** Displays current charges for the billing period as **\$0.00** and the estimated total for the billing period (Jun 1, 2020 - Jun 30, 2020) as **\$0.00**.
- Connections (2):** A red box highlights this section, which lists two connections: **node-red-kkhis-cloudant-1589315652695** and **Internet of Things Platform-uy-57641**.

11. Setup the Node RED editor.

The screenshot shows the Node-RED editor setup screen. The main content area contains the following text:


Welcome to your new Node-RED instance on IBM Cloud

We know you're eager to start wiring up your flows, but first there are a couple of tasks you should do:

- Secure your Node-RED editor
- Learn how to install additional nodes

At the bottom, there is a progress bar with four steps, the first of which is completed. Below the progress bar are **Previous** and **Next** buttons.

12. Continue clicking on Next and click on Finish at the end to launch the Node RED editor.



13. Click on Node-RED flow editor.

14. The Node-Red Editor window opens:

Before importing the flow following has to be installed:-

15. Install-


1. [node-red-contrib-ibm-watson-iot](#)
2. [node-red-contrib-scx-ibmiotapp](#)
3. [node-red-dashboard](#)
4. [node-red-node-twilio](#)
5. [node-red-contrib-ibmpush](#)

Incase it fails, add the following to package.json file

1. ["node-red-contrib-ibm-watson-iot": "0.x",](#)
2. ["node-red-contrib-scx-ibmiotapp": "0.x",](#)
3. ["node-red-dashboard": "2.x",](#)
4. ["node-red-node-twilio": "0.x",](#)
5. ["node-red-contrib-ibmpush": "0.x",](#)

Watson IoT Platform

Go to top menu ≡ -> Resource list->IoT Platform.



Launch the platform.

Press the + Add Device and create a Device id-"your_choice" and Device type-"esp8266". Skip over and press next until the token is prompted. Give one of your choice and you will finally reach the device drilldown page.

Note down the Org-ID, Device Type, Device ID and Authentication Token. It cannot be retrieved once you leave the page.

MQTT

We use MQTT to send data to IoT platform.

Add the following code in the Arduino IDE:

```
#define MQTT_HOST "ORG_ID.messaging.internetofthings.ibmcloud.com"
#define MQTT_PORT 1883
#define MQTT_DEVICEID "d:ORG_ID:DEVICE_TYPE:DEVICE_ID"
#define MQTT_USER "use-token-auth"
#define MQTT_TOKEN " //TOKEN THAT YOU ENTERED WHILE ADDING DEVICE IN IoT PLATFORM"
#define MQTT_TOPIC "iot-2/evt/status/fmt/json"
#define MQTT_TOPIC_DISPLAY "iot-2/cmd/display/fmt/json"
```

In place of ORG_ID, DEVIE_ID, DEVICE_TYPE and TOKEN, add the data you stored from the device drilldown page.

Open Connection Security Policy. Press the pencil icon next to Connection Security to edit the settings. Change the Default Security Level to TLS Optional, accept the Warning message by pressing the Ok button, then Save the change. This will use unsecured connection.

The rest of the mqtt code can be copied from the arduino code given.


On flashing the code, the data will be sent to the IoT platform.

Link to the code-” “

Generating your own API Key

Go to IoT platform. Go to the App section and click generate API Key.

The API Key is used as the username when connecting and the API Token is the password.



Setting up the nodes

Following nodes are installed in Node-RED

- **Node-red-contrib-scx-ibmiotapp and Node-red-contrib-ibm-watson iot**



- **Node-red-dashboard**

The following nodes are installed in the Node-RED to provide widgets that show up in your application user interface (UI).


▼ dashboard



- **Node-red-node-twilio**-twilio


Twilio APIs can add capabilities like voice, video, and messaging to applications.

▼ mobile




FLOW OF NODE-RED


- **Body Temperature**




- **pulse sensor flow**



- **Obstacle detection**



- **RFID tag**



IBM IoT input node



1. Drag an input > ibmiot in node from the palette on the left into the flow (the large blank area).
2. Double-click the new item.

Edit ibmiot in node

Delete Cancel Done

Properties

- Authentication: API Key
- API Key: IoT
- Input Type: Device Event
- Device Type: ☒ All or +
- Device Id: ☒ All or device id e.g. ab12cd231a21
- Event: ☒ All or +
- Format: ☒ All or json
- QoS: 0
- Name: IBM IoT
- Service: registered

Use the Input Type property to configure this node to receive Events sent by IoT Devices, Commands sent to IoT Devices, Status Messages referring to IoT Devices, or Status Messages referring to

☐ Enabled

3. Specify these properties:
 - Authentication: API Key
 - API Key: Add new ibmiot
 - Input Type: Device Event
 - Device Type: All
 - Device Id: All
 - Event: All
 - Format: All
4. Click the pencil icon.
5. Create an IoT with these parameters:
 - Name: "your_choice"
 - API Key and API Token: The values you created for the API Key previously
 - Server-Name: The host name you used (six characters for your specific ID followed by .messaging.internetofthings.ibmcloud.com)
6. Click Add and then Done

Change Node



Change node allow you to change a message payload or add new properties. You can use this node to affect the properties in a message, either by changing existing ones, deleting them or adding new properties.

1. Drag and drop and change node
2. Double click and add properties based on the input to the node

Template Node



1. Drag and drop template node
2. Double click and add the template to display if a value goes above or below a threshold.

Function Node



Function node are used to perform a specific function.

1. Drag and drop Function node
2. Double click on the node
3. Enter the code for the function to perform






Chart and gauge node



These node are used to display the obtained values in dashboard

1. Drag and drop the chart or gauge node.
2. Double click on the node
3. Put the group name is set as the name of the tab the graph or gauge belongs
4. Gauge node can be used to set max and min value for the input and can also be used to set the units of the input.

Debug Node




Debug nodes are used to check errors and values that are sent to different nodes.

1. Drag and drop debug node

2. Double click on the node
3. Error message if any will be shown in the debug messages

The 'Edit debug node' dialog box is shown. It has a title bar 'Edit debug node' and three buttons: 'Delete', 'Cancel', and 'Done'. Below the buttons is a 'Properties' section with a gear icon, a document icon, and a monitor icon. The 'Output' dropdown is set to 'complete msg object'. The 'To' section has two checkboxes: 'debug window' (checked) and 'system console' (unchecked). The 'Name' field is set to 'msg'.

Twilio out node



Add Account SID ,the number from which you will receive the call and Token. Also specify output as either sms or call and the number to which you will receive the alert.

The 'Edit twilio out node' dialog box is shown. It has a title bar 'Edit twilio out node' and three buttons: 'Delete', 'Cancel', and 'Done'. Below the buttons is a 'Properties' section with a gear icon, a document icon, and a monitor icon. The 'Twilio' field is set to 'Sms alert to admin'. The 'Output' dropdown is set to 'SMS'. The 'To' field is highlighted with a red box and contains a redacted phone number. The 'Name' field is set to 'Sms alert to admin'.

How to use CoviCare ?

The user and admin can visit website:

Link : <https://covicare.github.io/Repo/>


Demonstration video link :


https://drive.google.com/file/d/1jLFMDdcOZqKyhQXwOsE_FPvkyYnHIKWJ/view?usp=sharing

The user can either view the dashboard by logging in or can visit the Employment registration to register for the same.

For the purpose of attaining our goal to help as many people as we could, the Service is implemented as a user friendly website on a custom domain configured using a leading software development platform such as Github.

This ensures widening the reach of our service to millions of people through a widely supported and accessible hosting platform in a variety of browsers including Google Chrome, Mozilla Firefox and Microsoft edge .





We have also taken great care regarding the Friendliness and Accessibility of our website by testing it against the standards stated by the Wave Accessibility Evaluation tool [<https://wave.webaim.org>] and receiving acceptable results.


- Use the link provided to reach the welcome page of our website. Scroll down to find a brief description regarding our service.
- Click on the **Login** option . You will be taken to a new page with a login button in the top left corner.
- Click on this button to get a login form . Give the credentials: Sample username and password given below

Username :a

Password :a

The image shows a login form titled 'CoviCare Login' in the top left corner. A green 'Login' button is located in the top left of the form area. In the top right corner, there is a close button (an 'x' icon). The main content of the form is a large oval image of a wooden sign that says 'WELCOME'. Below this image, there are two input fields: 'Username' with a placeholder 'Enter Username' and 'Password' with a placeholder 'Enter Password'. Below the password field, there is a 'Login' button and a checkbox labeled 'Remember me'. At the bottom of the form, there is a red 'Cancel' button.

- Be careful with this form as u get only **three attempts** to login .On successful login,you will be taken to the node red dashboard page of Sevak.
- This page is complete with the visualisations and alert systems to monitor the health conditions and surroundings of the employee.
- On clicking the sandwich button,you can navigate to different tabs containing information about different parameters measured.



- The dashboard tabs can be selected and viewed from the list of tabs visible on clicking the option available in the top left corner.the list of tabs include:
 - Body Temperature [of the employee worker]
 - Pulse rate [pulse rate of the employee worker].
 - Distance monitoring
 - RFID tag
 - Current Covid Status
- All the above mentioned monitoring systems are equipped with an automatic alert system which
 - Location of the worker is monitored
 - Sends **sms or calls** the number provided with required information to the company in case symptoms arise in the worker.

The user can either view the dashboard by logging in or can visit the Employment registration to register for the same.

CLOUD LOGIN CREDENTIALS

Main cloud account:

Username: anjanaanil.mec@gmail.com

Password: Sevak@123

CONCLUSION

CoviCare is not just a tag. It is an all in one module that creates a hassle free work environment. From maintaining social distancing to monitoring your health, it will thoroughly ensure your well-being as well as those in your vicinity. We believe that CoviCare has the potential to become an essential part of everyday office life.